



ZYMO RESEARCH

The Beauty of Science is to Make Things Simple

04/14/2020

Kit Components

Product code	Description
--------------	-------------

D4010,D4011,D4064,D4065	Genomic DNA Clean & Concentrator
--------------------------------	---

Components:

D5201-1-50	ChIP DNA Binding Buffer
D4003-2-6	DNA Wash Buffer (Concentrate)
D3004-4-1	DNA Elution Buffer



Safety Data Sheet
acc. to OSHA HCS

Printing date 04/14/2020

Reviewed on 12/28/2015

1 Identification

- **Product identifier**
- **Trade name:** ChIP DNA Binding Buffer
- **Article number:** D5201-1-50, D5201-1-100
- **Application of the substance / the mixture** Laboratory Reagent
- **Details of the supplier of the safety data sheet**
- **Manufacturer/Supplier:**
Zymo Research Corp.
17062 Murphy Ave.
Irvine, CA 92614
USA
Phone: 1-949-679-1190 or 1-888-882-9682
sds@zymoresearch.com
- **Information department:** Product Safety Dept.
- **Emergency telephone number:**
During normal business hours (8 am to 5 pm Pacific Standard Time): +1 (949) 679 1190

2 Hazard(s) identification

- **Classification of the substance or mixture**



GHS07

- Acute Tox. 4 H302 Harmful if swallowed.
- Skin Irrit. 2 H315 Causes skin irritation.
- Eye Irrit. 2A H319 Causes serious eye irritation.

- **Label elements**
- **GHS label elements** The product is classified and labeled according to the Globally Harmonized System (GHS).
- **Hazard pictograms** GHS07
- **Signal word** Warning
- **Hazard-determining components of labeling:**
guanidinium chloride
- **Hazard statements**
Harmful if swallowed.
Causes skin irritation.
Causes serious eye irritation.
- **Precautionary statements**
Wash thoroughly after handling.
Do not eat, drink or smoke when using this product.
Wear protective gloves / eye protection / face protection.
If swallowed: Call a poison center/doctor if you feel unwell.
Rinse mouth.

(Contd. on page 2)



Safety Data Sheet acc. to OSHA HCS

Printing date 04/14/2020

Reviewed on 12/28/2015

Trade name: ChIP DNA Binding Buffer

(Contd. of page 1)

If on skin: Wash with plenty of water. Specific treatment (see on this label). If in eyes: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing. If skin irritation occurs: Get medical advice/attention. Take off contaminated clothing and wash it before reuse. If eye irritation persists: Get medical advice/attention. Dispose of contents/container in accordance with local/regional/national/international regulations.

Classification system:

NFPA ratings (scale 0 - 4)



Health = 2, Fire = 0, Reactivity = 0

HMIS-ratings (scale 0 - 4)



Health = 2, Fire = 0, Reactivity = 0

- Other hazards
Results of PBT and vPvB assessment
PBT: Not applicable.
vPvB: Not applicable.

3 Composition/information on ingredients

- Chemical characterization: Mixtures
Description: Mixture of the substances listed below with nonhazardous additions.

Dangerous components:

Table with 2 columns: CAS: 50-01-1, guanidinium chloride, and ≤30%

4 First-aid measures

- Description of first aid measures
General information: Symptoms of poisoning may even occur after several hours; therefore medical observation for at least 48 hours after the accident.
After inhalation: In case of unconsciousness place patient stably in side position for transportation.
After skin contact: Immediately wash with water and soap and rinse thoroughly.
After eye contact: Rinse opened eye for several minutes under running water. If symptoms persist, consult a doctor.
After swallowing: Do not induce vomiting; immediately call for medical help.

(Contd. on page 3)



Safety Data Sheet
acc. to OSHA HCS

Printing date 04/14/2020

Reviewed on 12/28/2015

Trade name: ChIP DNA Binding Buffer

(Contd. of page 2)

- **Information for doctor:**
- **Most important symptoms and effects, both acute and delayed** No further relevant information available.
- **Indication of any immediate medical attention and special treatment needed**
No further relevant information available.

5 Fire-fighting measures

- **Extinguishing media**
- **Suitable extinguishing agents:**
CO₂, extinguishing powder or water spray. Fight larger fires with water spray or alcohol resistant foam.
Use fire fighting measures that suit the environment.
- **Advice for firefighters**
- **Protective equipment:** Wear protective clothing.

6 Accidental release measures

- **Personal precautions, protective equipment and emergency procedures** Wear protective clothing.
- **Environmental precautions:** Dilute with plenty of water.
- **Methods and material for containment and cleaning up:**
Absorb with liquid-binding material (sand, diatomite, acid binders, universal binders, sawdust).
Dispose contaminated material as waste according to item 13.
Ensure adequate ventilation.
- **Reference to other sections**
See Section 7 for information on safe handling.
See Section 8 for information on personal protection equipment.
See Section 13 for disposal information.
- **Protective Action Criteria for Chemicals**

· **PAC-1:**
All components have the value 1.4 mg/m³.

· **PAC-2:**
All components have the value 16 mg/m³.

· **PAC-3:**
All components have the value 94 mg/m³.

7 Handling and storage

- **Handling:**
- **Precautions for safe handling**
Ensure good ventilation/exhaustion at the workplace.
Prevent formation of aerosols.

(Contd. on page 4)



Safety Data Sheet
acc. to OSHA HCS

Printing date 04/14/2020

Reviewed on 12/28/2015

Trade name: ChIP DNA Binding Buffer

(Contd. of page 3)

- Information about protection against explosions and fires: No special measures required.
Conditions for safe storage, including any incompatibilities
Storage:
Requirements to be met by storerooms and receptacles: No special requirements.
Information about storage in one common storage facility: Not required.
Further information about storage conditions: Keep receptacle tightly sealed.
Specific end use(s) Laboratory reagent

8 Exposure controls/personal protection

- Additional information about design of technical systems: No further data; see item 7.
Control parameters
Components with limit values that require monitoring at the workplace:
The product does not contain any relevant quantities of materials with critical values that have to be monitored at the workplace.
Additional information: The lists that were valid during the creation were used as basis.
Exposure controls
Personal protective equipment:
General protective and hygienic measures:
Keep away from foodstuffs, beverages and feed.
Immediately remove all soiled and contaminated clothing.
Wash hands before breaks and at the end of work.
Avoid contact with the eyes and skin.
Breathing equipment:
In case of brief exposure or low pollution use respiratory filter device. In case of intensive or longer exposure use respiratory protective device that is independent of circulating air.
Protection of hands:



Protective gloves

The glove material has to be impermeable and resistant to the product/ the substance/ the preparation. Due to missing tests no recommendation to the glove material can be given for the product/ the preparation/ the chemical mixture.

Selection of the glove material on consideration of the penetration times, rates of diffusion and the degradation

Material of gloves

The selection of the suitable gloves does not only depend on the material, but also on further marks of quality and varies from manufacturer to manufacturer. As the product is a preparation of several substances, the resistance of the glove material can not be calculated in advance and has therefore to be checked prior to the application.

Penetration time of glove material

The exact break through time has to be found out by the manufacturer of the protective gloves and has to be observed.

(Contd. on page 5)



Safety Data Sheet
acc. to OSHA HCS

Printing date 04/14/2020

Reviewed on 12/28/2015

Trade name: ChIP DNA Binding Buffer

(Contd. of page 4)

· **Eye protection:**



Tightly sealed goggles

9 Physical and chemical properties

· **Information on basic physical and chemical properties**

· **General Information**

· **Appearance:**

· Form:	Liquid
· Color:	Clear
· Odor:	Mild
· Odor threshold:	Not determined.

· **pH-value:** Not determined.

· **Change in condition**

· Melting point/Melting range:	Undetermined.
· Boiling point/Boiling range:	Undetermined.

· **Flash point:** Not applicable.

· **Flammability (solid, gaseous):** Not applicable.

· **Decomposition temperature:** Not determined.

· **Auto igniting:** Product is not selfigniting.

· **Danger of explosion:** Product does not present an explosion hazard.

· **Explosion limits:**

· Lower:	Not Applicable
· Upper:	Not Applicable

· **Vapor pressure:** Not determined.

· **Density:** Not determined.

· **Relative density:** Not determined.

· **Vapor density:** Not determined.

· **Evaporation rate:** Not determined.

· **Solubility in / Miscibility with**

· **Water:** Fully miscible.

· **Partition coefficient (n-octanol/water):** Not determined.

(Contd. on page 6)



Safety Data Sheet
acc. to OSHA HCS

Printing date 04/14/2020

Reviewed on 12/28/2015

Trade name: ChIP DNA Binding Buffer

(Contd. of page 5)

Table with 2 columns: Property (Viscosity, Solvent content, etc.) and Value (Not determined, 0.00%, etc.)

10 Stability and reactivity

- Reactivity, Chemical stability, Thermal decomposition, Possibility of hazardous reactions, Conditions to avoid, Incompatible materials, Hazardous decomposition products

11 Toxicological information

- Information on toxicological effects, Acute toxicity

LD/LC50 values that are relevant for classification:

CAS: 50-01-1 guanidinium chloride

Oral LD50 475 mg/kg (rat)

- Primary irritant effect: on the skin, on the eye; Sensitization; Additional toxicological information: Harmful, Irritant

Carcinogenic categories

IARC (International Agency for Research on Cancer)

None of the ingredients is listed.

(Contd. on page 7)



Safety Data Sheet
acc. to OSHA HCS

Printing date 04/14/2020

Reviewed on 12/28/2015

Trade name: ChIP DNA Binding Buffer

(Contd. of page 6)

NTP (National Toxicology Program)

None of the ingredients is listed.

OSHA-Ca (Occupational Safety & Health Administration)

None of the ingredients is listed.

12 Ecological information

- Toxicity
Aquatic toxicity: No further relevant information available.
Persistence and degradability No further relevant information available.
Behavior in environmental systems:
Bioaccumulative potential No further relevant information available.
Mobility in soil No further relevant information available.
Additional ecological information:
General notes: Not hazardous for water.
Results of PBT and vPvB assessment
PBT: Not applicable.
vPvB: Not applicable.
Other adverse effects No further relevant information available.

13 Disposal considerations

- Waste treatment methods
Recommendation:
Dispose of contents in accordance with local/regional/national, and international recommendations.
Uncleaned packagings:
Recommendation:
Dispose of container in accordance with local/regional/national and international recommendations.
Recommended cleansing agent: Water, if necessary with cleansing agents.

14 Transport information

- UN-Number
DOT, ADN, IMDG, IATA not regulated
UN proper shipping name
DOT, ADN, IMDG, IATA not regulated

(Contd. on page 8)



Safety Data Sheet acc. to OSHA HCS

Printing date 04/14/2020

Reviewed on 12/28/2015

Trade name: ChIP DNA Binding Buffer

(Contd. of page 7)

Table with 2 columns: Hazard/Regulation and Status. Rows include Transport hazard class(es), DOT, ADN, IMDG, IATA Class, Packing group, Environmental hazards, Special precautions for user, Transport in bulk according to Annex II of MARPOL 73/78 and the IBC Code, and UN Model Regulation.

15 Regulatory information

Safety, health and environmental regulations/legislation specific for the substance or mixture Sara

Section 355 (extremely hazardous substances): None of the ingredients is listed.

Section 313 (Specific toxic chemical listings): None of the ingredients is listed.

TSCA (Toxic Substances Control Act): All components have the value ACTIVE.

Hazardous Air Pollutants: None of the ingredients is listed.

Proposition 65

Chemicals known to cause cancer: None of the ingredients is listed.

Chemicals known to cause reproductive toxicity for females: None of the ingredients is listed.

Chemicals known to cause reproductive toxicity for males: None of the ingredients is listed.

Chemicals known to cause developmental toxicity: None of the ingredients is listed.

(Contd. on page 9)



Safety Data Sheet
acc. to OSHA HCS

Printing date 04/14/2020

Reviewed on 12/28/2015

Trade name: ChIP DNA Binding Buffer

(Contd. of page 8)

· **Carcinogenic categories**

· **EPA (Environmental Protection Agency)**

None of the ingredients is listed.

· **TLV (Threshold Limit Value established by ACGIH)**

None of the ingredients is listed.

· **NIOSH-Ca (National Institute for Occupational Safety and Health)**

None of the ingredients is listed.

· **GHS label elements** The product is classified and labeled according to the Globally Harmonized System (GHS).

· **Hazard pictograms** GHS07

· **Signal word** Warning

· **Hazard-determining components of labeling:**

guanidinium chloride

· **Hazard statements**

Harmful if swallowed.

Causes skin irritation.

Causes serious eye irritation.

· **Precautionary statements**

Wash thoroughly after handling.

Do not eat, drink or smoke when using this product.

Wear protective gloves / eye protection / face protection.

If swallowed: Call a poison center/doctor if you feel unwell.

Rinse mouth.

If on skin: Wash with plenty of water.

Specific treatment (see on this label).

If in eyes: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do.

Continue rinsing.

If skin irritation occurs: Get medical advice/attention.

Take off contaminated clothing and wash it before reuse.

If eye irritation persists: Get medical advice/attention.

Dispose of contents/container in accordance with local/regional/national/international regulations.

· **Chemical safety assessment:** A Chemical Safety Assessment has not been carried out.

16 Other information

This information is based on our present knowledge. However, this shall not constitute a guarantee for any specific product features and shall not establish a legally valid contractual relationship.

· **Department issuing SDS:**

Zymo Research Corp.

Safety Department

17062 Murphy Ave.

Irvine, CA 92614

USA

Phone: 1-949-679-1190 or 1-888-882-9682

(Contd. on page 10)



ZYMO RESEARCH

The Beauty of Science is to Make Things Simple

Page 10/10

Safety Data Sheet
acc. to OSHA HCS

Printing date 04/14/2020

Reviewed on 12/28/2015

Trade name: ChIP DNA Binding Buffer

(Contd. of page 9)

· **Contact:** sds@zymoresearch.com

· **Date of preparation / last revision** 04/14/2020 / -

· **Abbreviations and acronyms:**

ADR: Accord européen sur le transport des marchandises dangereuses par Route (European Agreement concerning the International Carriage of Dangerous Goods by Road)

IMDG: International Maritime Code for Dangerous Goods

DOT: US Department of Transportation

IATA: International Air Transport Association

ACGIH: American Conference of Governmental Industrial Hygienists

EINECS: European Inventory of Existing Commercial Chemical Substances

ELINCS: European List of Notified Chemical Substances

CAS: Chemical Abstracts Service (division of the American Chemical Society)

NFPA: National Fire Protection Association (USA)

HMIS: Hazardous Materials Identification System (USA)

VOC: Volatile Organic Compounds (USA, EU)

LC50: Lethal concentration, 50 percent

LD50: Lethal dose, 50 percent

PBT: Persistent, Bioaccumulative and Toxic

vPvB: very Persistent and very Bioaccumulative

NIOSH: National Institute for Occupational Safety

OSHA: Occupational Safety & Health

TLV: Threshold Limit Value

PEL: Permissible Exposure Limit

REL: Recommended Exposure Limit

Acute Tox. 4: Acute toxicity – Category 4

Skin Irrit. 2: Skin corrosion/irritation – Category 2

Eye Irrit. 2A: Serious eye damage/eye irritation – Category 2A

US



ZYMO RESEARCH

The Beauty of Science is to Make Things Simple

Safety Data Sheet
acc. to OSHA HCS

Printing date 04/14/2020

Reviewed on 12/16/2016

1 Identification

- **Product identifier**
- **Trade name:** DNA Wash Buffer (Concentrate)
- **Article number:** D4003-2-6, D4003-2-24, D4003-2-48
- **Application of the substance / the mixture** Laboratory Reagent
- **Details of the supplier of the safety data sheet**
- **Manufacturer/Supplier:**
Zymo Research Corp.
17062 Murphy Ave.
Irvine, CA 92614
USA
Phone: 1-949-679-1190 or 1-888-882-9682
sds@zymoresearch.com
- **Information department:** Product Safety Dept.
- **Emergency telephone number:**
During normal business hours (8 am to 5 pm Pacific Standard Time): +1 (949) 679 1190

2 Hazard(s) identification

- **Classification of the substance or mixture**
The product is not classified, according to the Globally Harmonized System (GHS).

- **Label elements**
- **GHS label elements** Void
- **Hazard pictograms** Void
- **Signal word** Void
- **Hazard statements** Void
- **Classification system:**
- **NFPA ratings (scale 0 - 4)**



- **HMIS-ratings (scale 0 - 4)**



- **Other hazards**
- **Results of PBT and vPvB assessment**
- **PBT:** Not applicable.
- **vPvB:** Not applicable.



Safety Data Sheet
acc. to OSHA HCS

Printing date 04/14/2020

Reviewed on 12/16/2016

Trade name: DNA Wash Buffer (Concentrate)

(Contd. of page 1)

3 Composition/information on ingredients

- **Chemical characterization: Mixtures**
- **Description:** Mixture of the substances listed below with nonhazardous additions.
- **Dangerous components:** Void

4 First-aid measures

- **Description of first aid measures**
- **General information:** No special measures required.
- **After inhalation:** Supply fresh air; consult doctor in case of complaints.
- **After skin contact:** Generally the product does not irritate the skin.
- **After eye contact:**
Rinse opened eye for several minutes under running water. Then consult a doctor.
Rinse opened eye for several minutes under running water.
- **After swallowing:** Do not induce vomiting; immediately call for medical help.
- **Information for doctor:**
- **Most important symptoms and effects, both acute and delayed** No further relevant information available.
- **Indication of any immediate medical attention and special treatment needed**
No further relevant information available.

5 Fire-fighting measures

- **Extinguishing media**
- **Suitable extinguishing agents:**
CO₂, extinguishing powder or water spray. Fight larger fires with water spray or alcohol resistant foam.
Use fire fighting measures that suit the environment.
- **Advice for firefighters**
- **Protective equipment:** Wear protective clothing.

6 Accidental release measures

- **Personal precautions, protective equipment and emergency procedures** Wear protective clothing.
- **Environmental precautions:** Dilute with plenty of water.
- **Methods and material for containment and cleaning up:**
Absorb with liquid-binding material (sand, diatomite, acid binders, universal binders, sawdust).
- **Reference to other sections**
See Section 7 for information on safe handling.
See Section 8 for information on personal protection equipment.
See Section 13 for disposal information.

(Contd. on page 3)



Safety Data Sheet
acc. to OSHA HCS

Printing date 04/14/2020

Reviewed on 12/16/2016

Trade name: DNA Wash Buffer (Concentrate)

(Contd. of page 2)

Protective Action Criteria for Chemicals

PAC-1:

Table with 3 columns: CAS number, chemical name, and concentration (mg/m³). Rows include 2-amino-2-(hydroxymethyl)propane-1,3-diolhydrochloride and Edetate Disodium, Dihydrate.

PAC-2:

Table with 3 columns: CAS number, chemical name, and concentration (mg/m³). Rows include 2-amino-2-(hydroxymethyl)propane-1,3-diolhydrochloride and Edetate Disodium, Dihydrate.

PAC-3:

Table with 3 columns: CAS number, chemical name, and concentration (mg/m³). Rows include 2-amino-2-(hydroxymethyl)propane-1,3-diolhydrochloride and Edetate Disodium, Dihydrate.

7 Handling and storage

- Handling:
Precautions for safe handling: No special measures required.
Information about protection against explosions and fires: No special measures required.
Conditions for safe storage, including any incompatibilities
Storage:
Requirements to be met by storerooms and receptacles: No special requirements.
Information about storage in one common storage facility: Not required.
Further information about storage conditions: None.
Specific end use(s): Laboratory reagent

8 Exposure controls/personal protection

- Additional information about design of technical systems: No further data; see item 7.
Control parameters
Components with limit values that require monitoring at the workplace:
The product does not contain any relevant quantities of materials with critical values that have to be monitored at the workplace.
Additional information: The lists that were valid during the creation were used as basis.
Exposure controls
Personal protective equipment:
General protective and hygienic measures:
The usual precautionary measures for handling chemicals should be followed.
Breathing equipment: Not required.
Protection of hands:
The glove material has to be impermeable and resistant to the product/ the substance/ the preparation.

(Contd. on page 4)



Safety Data Sheet
acc. to OSHA HCS

Printing date 04/14/2020

Reviewed on 12/16/2016

Trade name: DNA Wash Buffer (Concentrate)

(Contd. of page 3)

Due to missing tests no recommendation to the glove material can be given for the product/ the preparation/ the chemical mixture.

Selection of the glove material on consideration of the penetration times, rates of diffusion and the degradation

Material of gloves

The selection of the suitable gloves does not only depend on the material, but also on further marks of quality and varies from manufacturer to manufacturer. As the product is a preparation of several substances, the resistance of the glove material can not be calculated in advance and has therefore to be checked prior to the application.

Penetration time of glove material

The exact break through time has to be found out by the manufacturer of the protective gloves and has to be observed.

Eye protection: Goggles recommended during refilling.

9 Physical and chemical properties

Information on basic physical and chemical properties

General Information

Appearance:

Table with 2 columns: Property (Form, Color, Odor, Odor threshold) and Value (Liquid, Clear, Odorless, Not determined).

pH-value: Not determined.

Change in condition

Table with 2 columns: Property (Melting point/Melting range, Boiling point/Boiling range) and Value (Undetermined, Undetermined).

Flash point: Not applicable.

Flammability (solid, gaseous): Not applicable.

Decomposition temperature: Not determined.

Auto igniting: Product is not selfigniting.

Danger of explosion: Product does not present an explosion hazard.

Explosion limits:

Table with 2 columns: Property (Lower, Upper) and Value (Not Applicable, Not Applicable).

Vapor pressure: Not determined.

Table with 2 columns: Property (Density, Relative density, Vapor density, Evaporation rate) and Value (Not determined, Not determined, Not determined, Not determined).

(Contd. on page 5)



Safety Data Sheet
acc. to OSHA HCS

Printing date 04/14/2020

Reviewed on 12/16/2016

Trade name: DNA Wash Buffer (Concentrate)

(Contd. of page 4)

Table with 2 columns: Property and Value. Rows include Solubility in / Miscibility with Water, Partition coefficient, Viscosity, Solvent content, and Solids content.

10 Stability and reactivity

- Reactivity, Chemical stability, Thermal decomposition / conditions to be avoided, Possibility of hazardous reactions, Conditions to avoid, Incompatible materials, Hazardous decomposition products.

11 Toxicological information

- Information on toxicological effects, Acute toxicity, Primary irritant effect, Sensitization, Additional toxicological information, Carcinogenic categories.

Table with 1 column: IARC (International Agency for Research on Cancer). Value: None of the ingredients is listed.

(Contd. on page 6)



Safety Data Sheet
acc. to OSHA HCS

Printing date 04/14/2020

Reviewed on 12/16/2016

Trade name: DNA Wash Buffer (Concentrate)

(Contd. of page 5)

NTP (National Toxicology Program)

None of the ingredients is listed.

OSHA-Ca (Occupational Safety & Health Administration)

None of the ingredients is listed.

12 Ecological information

- Toxicity
Aquatic toxicity: No further relevant information available.
Persistence and degradability No further relevant information available.
Behavior in environmental systems:
Bioaccumulative potential No further relevant information available.
Mobility in soil No further relevant information available.
Additional ecological information:
General notes: Not hazardous for water.
Results of PBT and vPvB assessment
PBT: Not applicable.
vPvB: Not applicable.
Other adverse effects No further relevant information available.

13 Disposal considerations

- Waste treatment methods
Recommendation: Smaller quantities can be disposed of with household waste.
Uncleaned packagings:
Recommendation:
Dispose of container in accordance with local/regional/national and international recommendations.
Recommended cleansing agent: Water, if necessary with cleansing agents.

14 Transport information

- UN-Number
DOT, ADN, IMDG, IATA not regulated
UN proper shipping name
DOT, ADN, IMDG, IATA not regulated
Transport hazard class(es)
DOT, ADN, IMDG, IATA
Class not regulated

(Contd. on page 7)



Safety Data Sheet
acc. to OSHA HCS

Printing date 04/14/2020

Reviewed on 12/16/2016

Trade name: DNA Wash Buffer (Concentrate)

(Contd. of page 6)

Table with 2 columns: Hazard/Regulation and Status. Rows include Packing group, DOT, IMDG, IATA, Environmental hazards, Special precautions for user, Transport in bulk according to Annex II of MARPOL73/78 and the IBC Code, and UN Model Regulation.

15 Regulatory information

Safety, health and environmental regulations/legislation specific for the substance or mixture
Sara

Table with 2 columns: Regulation/Category and Status. Rows include Section 355, Section 313, TSCA, Hazardous Air Pollutants, Proposition 65, Chemicals known to cause cancer, Chemicals known to cause reproductive toxicity for females, Chemicals known to cause reproductive toxicity for males, Chemicals known to cause developmental toxicity, Carcinogenic categories, EPA, and TLV.

(Contd. on page 8)



ZYMO RESEARCH

The Beauty of Science is to Make Things Simple

Page 8/8

Safety Data Sheet *acc. to OSHA HCS*

Printing date 04/14/2020

Reviewed on 12/16/2016

Trade name: DNA Wash Buffer (Concentrate)

(Contd. of page 7)

· **NIOSH-Ca (National Institute for Occupational Safety and Health)**

None of the ingredients is listed.

· **GHS label elements** Void

· **Hazard pictograms** Void

· **Signal word** Void

· **Hazard statements** Void

· **Chemical safety assessment:** A Chemical Safety Assessment has not been carried out.

16 Other information

This information is based on our present knowledge. However, this shall not constitute a guarantee for any specific product features and shall not establish a legally valid contractual relationship.

· **Department issuing SDS:**

Zymo Research Corp.

Safety Department

17062 Murphy Ave.

Irvine, CA 92614

USA

Phone: 1-949-679-1190 or 1-888-882-9682

· **Contact:** sds@zymoresearch.com

· **Date of preparation / last revision** 04/14/2020 / -

· **Abbreviations and acronyms:**

ADR: Accord européen sur le transport des marchandises dangereuses par Route (European Agreement concerning the International Carriage of Dangerous Goods by Road)

IMDG: International Maritime Code for Dangerous Goods

DOT: US Department of Transportation

IATA: International Air Transport Association

ACGIH: American Conference of Governmental Industrial Hygienists

EINECS: European Inventory of Existing Commercial Chemical Substances

ELINCS: European List of Notified Chemical Substances

CAS: Chemical Abstracts Service (division of the American Chemical Society)

NFPA: National Fire Protection Association (USA)

HMIS: Hazardous Materials Identification System (USA)

VOC: Volatile Organic Compounds (USA, EU)

PBT: Persistent, Bioaccumulative and Toxic

vPvB: very Persistent and very Bioaccumulative

NIOSH: National Institute for Occupational Safety

OSHA: Occupational Safety & Health

TLV: Threshold Limit Value

PEL: Permissible Exposure Limit

REL: Recommended Exposure Limit

US



Safety Data Sheet
acc. to OSHA HCS

Printing date 04/14/2020

Reviewed on 12/03/2015

1 Identification

- **Product identifier**
- **Trade name:** DNA Elution Buffer
- **Article number:** D3004-4-1, D3004-4-4, D3004-4-10, D3004-4-16, D3004-4-50
- **Application of the substance / the mixture** Laboratory Reagent
- **Details of the supplier of the safety data sheet**
- **Manufacturer/Supplier:**
Zymo Research Corp.
17062 Murphy Ave.
Irvine, CA 92614
USA
Phone: 1-949-679-1190 or 1-888-882-9682
sds@zymoresearch.com
- **Information department:** Product Safety Dept.
- **Emergency telephone number:**
During normal business hours (8 am to 5 pm Pacific Standard Time): +1 (949) 679 1190

2 Hazard(s) identification

- **Classification of the substance or mixture**
The product is not classified, according to the Globally Harmonized System (GHS).

- **Label elements**
- **GHS label elements** Void
- **Hazard pictograms** Void
- **Signal word** Void
- **Hazard statements** Void
- **Classification system:**
- **NFPA ratings (scale 0 - 4)**



- **HMIS-ratings (scale 0 - 4)**



- **Other hazards**
- **Results of PBT and vPvB assessment**
- **PBT:** Not applicable.
- **vPvB:** Not applicable.



Safety Data Sheet
acc. to OSHA HCS

Printing date 04/14/2020

Reviewed on 12/03/2015

Trade name: DNA Elution Buffer

(Contd. of page 1)

3 Composition/information on ingredients

- **Chemical characterization: Mixtures**
- **Description:** Mixture of the substances listed below with nonhazardous additions.
- **Dangerous components:** Void

4 First-aid measures

- **Description of first aid measures**
- **General information:** No special measures required.
- **After inhalation:** Supply fresh air; consult doctor in case of complaints.
- **After skin contact:** Generally the product does not irritate the skin.
- **After eye contact:**
Rinse opened eye for several minutes under running water. Then consult a doctor.
Rinse opened eye for several minutes under running water.
- **After swallowing:** Do not induce vomiting; immediately call for medical help.
- **Information for doctor:**
- **Most important symptoms and effects, both acute and delayed** No further relevant information available.
- **Indication of any immediate medical attention and special treatment needed**
No further relevant information available.

5 Fire-fighting measures

- **Extinguishing media**
- **Suitable extinguishing agents:**
CO₂, extinguishing powder or water spray. Fight larger fires with water spray or alcohol resistant foam.
Use fire fighting measures that suit the environment.
- **Advice for firefighters**
- **Protective equipment:** Wear protective clothing.

6 Accidental release measures

- **Personal precautions, protective equipment and emergency procedures** Wear protective clothing.
- **Environmental precautions:**
Dilute with plenty of water.
Do not allow to enter sewers/ surface or ground water.
- **Methods and material for containment and cleaning up:**
Absorb with liquid-binding material (sand, diatomite, acid binders, universal binders, sawdust).
- **Reference to other sections**
See Section 7 for information on safe handling.
See Section 8 for information on personal protection equipment.
See Section 13 for disposal information.

(Contd. on page 3)



Safety Data Sheet
acc. to OSHA HCS

Printing date 04/14/2020

Reviewed on 12/03/2015

Trade name: DNA Elution Buffer

(Contd. of page 2)

Protective Action Criteria for Chemicals

PAC-1:

Table with 3 columns: CAS number, chemical name, and concentration. Rows include CAS: 1185-53-1 for 2-amino-2-(hydroxymethyl)propane-1,3-diolhydrochloride at 12 mg/m³ and CAS: 6381-92-6 for Edetate Disodium, Dihydrate at 30 mg/m³.

PAC-2:

Table with 3 columns: CAS number, chemical name, and concentration. Rows include CAS: 1185-53-1 for 2-amino-2-(hydroxymethyl)propane-1,3-diolhydrochloride at 130 mg/m³ and CAS: 6381-92-6 for Edetate Disodium, Dihydrate at 330 mg/m³.

PAC-3:

Table with 3 columns: CAS number, chemical name, and concentration. Rows include CAS: 1185-53-1 for 2-amino-2-(hydroxymethyl)propane-1,3-diolhydrochloride at 790 mg/m³ and CAS: 6381-92-6 for Edetate Disodium, Dihydrate at 2,000 mg/m³.

7 Handling and storage

- Handling:
Precautions for safe handling: No special measures required.
Information about protection against explosions and fires: No special measures required.
Conditions for safe storage, including any incompatibilities
Storage:
Requirements to be met by storerooms and receptacles: No special requirements.
Information about storage in one common storage facility: Not required.
Further information about storage conditions: None.
Specific end use(s): Laboratory reagent

8 Exposure controls/personal protection

- Additional information about design of technical systems: No further data; see item 7.
Control parameters
Components with limit values that require monitoring at the workplace:
The product does not contain any relevant quantities of materials with critical values that have to be monitored at the workplace.
Additional information: The lists that were valid during the creation were used as basis.
Exposure controls
Personal protective equipment:
General protective and hygienic measures:
The usual precautionary measures for handling chemicals should be followed.
Breathing equipment: Not required.
Protection of hands:
The glove material has to be impermeable and resistant to the product/ the substance/ the preparation.

(Contd. on page 4)



Safety Data Sheet
acc. to OSHA HCS

Printing date 04/14/2020

Reviewed on 12/03/2015

Trade name: DNA Elution Buffer

(Contd. of page 3)

Due to missing tests no recommendation to the glove material can be given for the product/ the preparation/ the chemical mixture.

Selection of the glove material on consideration of the penetration times, rates of diffusion and the degradation

Material of gloves

The selection of the suitable gloves does not only depend on the material, but also on further marks of quality and varies from manufacturer to manufacturer. As the product is a preparation of several substances, the resistance of the glove material can not be calculated in advance and has therefore to be checked prior to the application.

Penetration time of glove material

The exact break through time has to be found out by the manufacturer of the protective gloves and has to be observed.

Eye protection: Goggles recommended during refilling.

9 Physical and chemical properties

Information on basic physical and chemical properties

General Information

Appearance:

Table with 2 columns: Property (Form, Color, Odor, Odor threshold) and Value (Liquid, Clear, Odorless, Not determined).

pH-value: Not determined.

Change in condition

Table with 2 columns: Property (Melting point/Melting range, Boiling point/Boiling range) and Value (Undetermined, Undetermined).

Flash point: Not applicable.

Flammability (solid, gaseous): Not applicable.

Decomposition temperature: Not determined.

Auto igniting: Product is not selfigniting.

Danger of explosion: Product does not present an explosion hazard.

Explosion limits:

Table with 2 columns: Property (Lower, Upper) and Value (Not Applicable, Not Applicable).

Vapor pressure: Not determined.

Table with 2 columns: Property (Density, Relative density, Vapor density, Evaporation rate) and Value (Not determined, Not determined, Not determined, Not determined).

(Contd. on page 5)



Safety Data Sheet
acc. to OSHA HCS

Printing date 04/14/2020

Reviewed on 12/03/2015

Trade name: DNA Elution Buffer

(Contd. of page 4)

· Solubility in / Miscibility with Water:	Fully miscible.
· Partition coefficient (n-octanol/water):	Not determined.
· Viscosity:	
· Dynamic:	Not determined.
· Kinematic:	Not determined.
· Solvent content:	
· VOC content:	0.00 % 0.0 g/l / 0.00 lb/gal
· Solids content:	2.0 %
· Other information	No further relevant information available.

10 Stability and reactivity

- **Reactivity** No further relevant information available.
- **Chemical stability**
- **Thermal decomposition / conditions to be avoided:** No decomposition if used according to specifications.
- **Possibility of hazardous reactions** No dangerous reactions known.
- **Conditions to avoid** No further relevant information available.
- **Incompatible materials:** No further relevant information available.
- **Hazardous decomposition products:** No dangerous decomposition products known.

11 Toxicological information

- **Information on toxicological effects**
- **Acute toxicity:**
- **Primary irritant effect:**
- **on the skin:** No irritant effect.
- **on the eye:** No irritating effect.
- **Sensitization:** No sensitizing effects known.
- **Additional toxicological information:**
The product is not subject to classification according to internally approved calculation methods for preparations: When used and handled according to specifications, the product does not have any harmful effects according to our experience and the information provided to us.
- **Carcinogenic categories**

· IARC (International Agency for Research on Cancer)
None of the ingredients is listed.

(Contd. on page 6)



Safety Data Sheet acc. to OSHA HCS

Printing date 04/14/2020

Reviewed on 12/03/2015

Trade name: DNA Elution Buffer

(Contd. of page 5)

· NTP (National Toxicology Program)

None of the ingredients is listed.

· OSHA-Ca (Occupational Safety & Health Administration)

None of the ingredients is listed.

12 Ecological information

- Toxicity
· Aquatic toxicity: No further relevant information available.
· Persistence and degradability No further relevant information available.
· Behavior in environmental systems:
· Bioaccumulative potential No further relevant information available.
· Mobility in soil No further relevant information available.
· Additional ecological information:
· General notes:
Water hazard class 1 (Self-assessment): slightly hazardous for water
Do not allow undiluted product or large quantities of it to reach ground water, water course or sewage system.
· Results of PBT and vPvB assessment
· PBT: Not applicable.
· vPvB: Not applicable.
· Other adverse effects No further relevant information available.

13 Disposal considerations

- Waste treatment methods
· Recommendation: Smaller quantities can be disposed of with household waste.
· Uncleaned packagings:
· Recommendation:
Dispose of container in accordance with local/regional/national and international recommendations.
· Recommended cleansing agent: Water, if necessary with cleansing agents.

14 Transport information

- UN-Number
· DOT, ADN, IMDG, IATA not regulated
· UN proper shipping name
· DOT, ADN, IMDG, IATA not regulated

(Contd. on page 7)



Safety Data Sheet
acc. to OSHA HCS

Printing date 04/14/2020

Reviewed on 12/03/2015

Trade name: DNA Elution Buffer

(Contd. of page 6)

· Transport hazard class(es)	
· DOT, ADN, IMDG, IATA	
· Class	not regulated
· Packing group	
· DOT, IMDG, IATA	not regulated
· Environmental hazards:	Not applicable.
· Special precautions for user	Not applicable.
· Transport in bulk according to Annex II of MARPOL 73/78 and the IBC Code	Not applicable.
· UN "Model Regulation":	not regulated

15 Regulatory information

· **Safety, health and environmental regulations/legislation specific for the substance or mixture**
· **Sara**

· **Section 355 (extremely hazardous substances):**
None of the ingredients is listed.

· **Section 313 (Specific toxic chemical listings):**
None of the ingredients is listed.

· **TSCA (Toxic Substances Control Act):**
CAS: 1185-53-1 | 2-amino-2-(hydroxymethyl)propane-1,3-diolhydrochloride | ACTIVE

· **Hazardous Air Pollutants**
None of the ingredients is listed.

· **Proposition 65**

· **Chemicals known to cause cancer:**
None of the ingredients is listed.

· **Chemicals known to cause reproductive toxicity for females:**
None of the ingredients is listed.

· **Chemicals known to cause reproductive toxicity for males:**
None of the ingredients is listed.

· **Chemicals known to cause developmental toxicity:**
None of the ingredients is listed.

(Contd. on page 8)



Safety Data Sheet
acc. to OSHA HCS

Printing date 04/14/2020

Reviewed on 12/03/2015

Trade name: DNA Elution Buffer

(Contd. of page 7)

Carcinogenic categories

EPA (Environmental Protection Agency)

None of the ingredients is listed.

TLV (Threshold Limit Value established by ACGIH)

None of the ingredients is listed.

NIOSH-Ca (National Institute for Occupational Safety and Health)

None of the ingredients is listed.

GHS label elements Void

Hazard pictograms Void

Signal word Void

Hazard statements Void

Chemical safety assessment: A Chemical Safety Assessment has not been carried out.

16 Other information

This information is based on our present knowledge. However, this shall not constitute a guarantee for any specific product features and shall not establish a legally valid contractual relationship.

Department issuing SDS:

Zymo Research Corp.
Safety Department
17062 Murphy Ave.
Irvine, CA 92614
USA

Phone: 1-949-679-1190 or 1-888-882-9682

Contact: sds@zymoresearch.com

Date of preparation / last revision 04/14/2020 / -

Abbreviations and acronyms:

- ADR: Accord européen sur le transport des marchandises dangereuses par Route (European Agreement concerning the International Carriage of Dangerous Goods by Road)
IMDG: International Maritime Code for Dangerous Goods
DOT: US Department of Transportation
IATA: International Air Transport Association
ACGIH: American Conference of Governmental Industrial Hygienists
EINECS: European Inventory of Existing Commercial Chemical Substances
ELINCS: European List of Notified Chemical Substances
CAS: Chemical Abstracts Service (division of the American Chemical Society)
NFPA: National Fire Protection Association (USA)
HMIS: Hazardous Materials Identification System (USA)
VOC: Volatile Organic Compounds (USA, EU)
PBT: Persistent, Bioaccumulative and Toxic
vPvB: very Persistent and very Bioaccumulative
NIOSH: National Institute for Occupational Safety
OSHA: Occupational Safety & Health
TLV: Threshold Limit Value
PEL: Permissible Exposure Limit
REL: Recommended Exposure Limit