



PROTECTION

# *Microflex<sup>®</sup> Gloves*

*Passionate for Science.*



[www.starlabgroup.com](http://www.starlabgroup.com)

# Microflex® Gloves

## Glove overview

STARLAB is committed to supplying the best hand protection to its customers. To provide a wider choice for our customers, we also sell a range of high quality Microflex® gloves for use in the laboratory to help you select the best glove for the job in hand.



**XCEED®**  
Ergonomic  
Blue Nitrile

**SUPRENO®**  
Blue  
Nitrile

**93-853**  
Extended  
Blue Nitrile

**93-243**  
Extra Extended  
Blue Nitrile

**73-847**  
Ergonomic  
Neoprene

**93-260**  
Chemical Resistant  
Triple Layered


GLOVE	Material	Interior	Exterior	Colour	Length	Thickness (palm)	AQL	PPE Cat III	High Risk	Pack Size
<b>XCEED®</b>	Nitrile	Powder-free	Textured Fingers	Blue	Standard 240 mm	0.070mm	0.65	✓	✗	10 x 250'
<b>SUPRENO®</b>	Nitrile	Powder-free	Textured Fingers	Blue	Standard 245 mm	0.110mm	0.65	✓	✗	10 x 100
<b>93-853</b>	Nitrile	Powder-free	Textured Fingers	Blue	Extended 295 mm	0.140mm	0.65	✓	✗	10 x 50
<b>93-243</b>	Nitrile	Powder-free	Textured Palm and Fingers	Blue	Extra extended 400 mm	0.120mm	1.5	✓	✗	10 x 100
<b>73-847</b>	Neoprene	Powder-free	Textured Fingers	Dark Green	Standard 245 mm	0.100mm	1.5	✓	✗	10 x 100
<b>93-260</b>	Nitrile and Neoprene	Powder-free	Textured Fingers	Green outside, Blue inside	Extended 300 mm	0.198mm	0.65	✓	✓	10 x 50

## Product features

- ▶ ASTM D6978 Approved for handling cytostatic agents (Supreno® & 93-853)
- ▶ EN 1149 (Part 1, 2 & 3) Antistatic behaviour (93-853 & 73-847)
- ▶ EN 388 Protection against mechanical risks (93-260)
- ▶ EN 421 Protective gloves against ionising radiation and radioactive contamination (93-260 & 73-847)
- ▶ Certified ergonomics (XCEED® & 73-847)
- ▶ Forearm protection/ Extended glove length > 295 mm (93-853,93-243,93-260)

**INFO:** For a complete list of standards applicable to each glove, please refer to the relevant product information.





EN ISO 374-1/Type B



EN ISO 374-5:2016



AQL 0.65

### XCEED® – Comfortable ergonomic nitrile protection

A combination of exceptional comfort and strength.

Recommended for use when:

- ▶ Work demands great dexterity and you need a glove to reduce hand fatigue
- ▶ You are looking for significantly stronger protection
- ▶ Working in environments that demand consistent grip and reliability
- ▶ Using Ethidium Bromide (Breakthrough: >480 mins)




**PPE Cat. III Complex Design**  
**Blue Nitrile Glove. Textured Fingers**  
**Length: 240 mm. 250 Gloves per Box<sup>1)</sup>**

ORDER INFORMATION			
Description	Pack Size	Cat. No.	Price
XCEED® Nitrile Gloves, XS	10 × 250	XC-INT-XS	£358.00
XCEED® Nitrile Gloves, S	10 × 250	XC-INT-S	£358.00
XCEED® Nitrile Gloves, M	10 × 250	XC-INT-M	£358.00
XCEED® Nitrile Gloves, L	10 × 250	XC-INT-L	£358.00
XCEED® Nitrile Gloves, XL <sup>1)</sup>	10 × 230	XC-INT-XL	£358.00

<sup>1)</sup> XL 230 gloves

[www.starlab.click/xceed](http://www.starlab.click/xceed)



EN ISO 374-1/Type B



EN ISO 374-5:2016



AQL 0.65

### SUPRENO® – Everyday nitrile protection

Comfortable nitrile protection for everyday use.

Recommended for use when:

- ▶ Contact with chemicals is a standard activity each day - SUPRENO® provides excellent chemical splash protection
- ▶ Strength and protection are of great importance
- ▶ Working in environments that demand consistent grip and reliability



**PPE Cat. III Complex Design**  
**Blue Nitrile Glove. Textured Fingers**  
**Length: 245 mm. 100 Gloves per Box**

ORDER INFORMATION			
Description	Pack Size	Cat. No.	Price
SUPRENO®, Nitrile Gloves, XS	10 × 100	SU-INT-XS	£155.35
SUPRENO®, Nitrile Gloves, S	10 × 100	SU-INT-S	£155.35
SUPRENO®, Nitrile Gloves, M	10 × 100	SU-INT-M	£155.35
SUPRENO®, Nitrile Gloves, L	10 × 100	SU-INT-L	£155.35
SUPRENO®, Nitrile Gloves, XL	10 × 100	SU-INT-XL	£155.35

[www.starlab.click/supreno](http://www.starlab.click/supreno)

# Microflex® Gloves



**AQL 0.65**

ORDER INFORMATION			
Description	Pack Size	Cat. No.	Price
Microflex® 93-853 Nitrile Glove, Ext. Cuff, XS	10 x 50	93-853-XS	£137.79
Microflex® 93-853 Nitrile Glove, Ext. Cuff, S	10 x 50	93-853-S	£137.79
Microflex® 93-853 Nitrile Glove, Ext. Cuff, M	10 x 50	93-853-M	£137.79
Microflex® 93-853 Nitrile Glove, Ext. Cuff, L	10 x 50	93-853-L	£137.79
Microflex® 93-853 Nitrile Glove, Ext. Cuff, XL	10 x 50	93-853-XL	£137.79

[www.starlab.click/93-853](http://www.starlab.click/93-853)

## 93-853 – Reliable hazard protection

Extended cuff provides additional protection over the wrist for hazardous environments. Advanced nitrile formulation offers up to three times the puncture resistance of latex or vinyl disposable gloves.

- ▶ Excellent chemical breakthrough resistance
- ▶ Over wrist protection against laboratory hazards; glove length 295 mm
- ▶ Very high barrier properties (AQL 0.65)
- ▶ Anti-static behavior
- ▶ Approved for use with chemotherapy drugs
- ▶ Medical examination glove

### PPE Cat. III Complex Design

**Blue Nitrile Glove. Textured Fingers**

**Length: 295 mm. 50 Gloves per Bag**



**AQL 1.5**

ORDER INFORMATION			
Description	Pack Size	Cat. No.	Price
Microflex® 93-243, Extra Long, XS	10 x 100	93-243-XS	£206.00
Microflex® 93-243, Extra Long, S	10 x 100	93-243-S	£206.00
Microflex® 93-243, Extra Long, M	10 x 100	93-243-M	£206.00
Microflex® 93-243, Extra Long, L	10 x 100	93-243-L	£206.00
Microflex® 93-243, Extra Long, XL	10 x 100	93-243-XL	£206.00
Microflex® 93-243, Extra Long, XXL	10 x 100	93-243-XXL	£206.00

[www.starlab.click/93-243](http://www.starlab.click/93-243)

## 93-243 – Superior forearm protection

Soft durable nitrile material ensures extended protection up to the elbow with a comfortable feel.

- ▶ Provides extra security over the wrist and forearm
- ▶ Unique cuff design prevents roll-down
- ▶ Textured palm and fingers for a confident secure grip




### PPE Cat. III Complex Design

**Blue Nitrile Glove. Textured Fingers**

**Length: 400 mm. 100 Gloves per Bag**



400 mm glove length for superior forearm protection



### 73-847 – Ergonomically certified

Soft and pliable neoprene formulation that provides effective resistance to many chemicals, and provides a comfortable and tactile wear that is close in feel to latex.


- ▶ Ergonomically certified for excellent comfort during use
- ▶ Excellent grip in both wet or dry conditions
- ▶ Textured fingers
- ▶ Anti-static behavior



**PPE Cat. III Complex Design**  
**Green Neoprene Glove. Textured Fingers**  
**Length: 245 mm. 100 Gloves per Box**

ORDER INFORMATION			
Description	Pack Size	Cat. No.	Price
Microflex® 73-847, Neoprene Glove, Ergonomic, XS	10 × 100	73-847-XS	£118.45
Microflex® 73-847, Neoprene Glove, Ergonomic, S	10 × 100	73-847-S	£118.45
Microflex® 73-847, Neoprene Glove, Ergonomic, M	10 × 100	73-847-M	£118.45
Microflex® 73-847, Neoprene Glove, Ergonomic, L	10 × 100	73-847-L	£118.45
Microflex® 73-847, Neoprene Glove, Ergonomic, XL	10 × 100	73-847-XL	£118.45

[www.starlab.click/73-847](http://www.starlab.click/73-847)



### 93-260 – Chemical protection and comfort

Innovative three layered glove design comprising nitrile and neoprene that provides superior protection against harsh chemicals, whilst still providing comfort and dexterity for the user.

- ▶ Ergonomic design for excellent comfort during use
- ▶ Lower acceptable pinhole rate (0.65 AQL) for reliable protection against hazardous substances
- ▶ Longer length for even greater protection
- ▶ Type A classification




**PPE Cat. III Complex Design**  
**Green Nitrile/Neoprene Glove. Textured Fingers**  
**Length: 300 mm. 50 Gloves per Box**

ORDER INFORMATION			
Description	Pack Size	Cat. No.	Price
Microflex® 93-260 Nitrile/Neoprene Glove, Chemical Resistant, XS	10 × 50	93-260-XS	£144.20
Microflex® 93-260 Nitrile/Neoprene Glove, Chemical Resistant, S	10 × 50	93-260-S	£144.20
Microflex® 93-260, Nitrile/Neoprene Glove, Chemical Resistant, M	10 × 50	93-260-M	£144.20
Microflex® 93-260 Nitrile/Neoprene Glove, Chemical Resistant, L	10 × 50	93-260-L	£144.20
Microflex® 93-260 Nitrile/Neoprene Glove, Chemical Resistant, XL	10 × 50	93-260-XL	£144.20
Microflex® 93-260 Nitrile/Neoprene Glove, Chemical Resistant, XXL	10 × 50	93-260-XXL	£144.20

[www.starlab.click/93-260](http://www.starlab.click/93-260)

Innovative three layer design





[www.starlab.co.uk](http://www.starlab.co.uk)

**STARLAB (UK), Ltd**

5 Tanners Drive  
Milton Keynes MK14 5BU  
United Kingdom

T: +44 (0)1908 283800  
F: +44 (0)1908 283802  
[infoline@starlab.co.uk](mailto:infoline@starlab.co.uk)  
[www.starlab.co.uk](http://www.starlab.co.uk)

**STARLAB International GmbH**

Neuer Höltigbaum 38  
22143 Hamburg  
Germany

T: +49 (0)40 675 99 39 0  
F: +49 (0)40 675 99 39 20  
[info@starlab.de](mailto:info@starlab.de)  
[www.starlabgroup.com](http://www.starlabgroup.com)

**STARLAB GmbH**

Neuer Höltigbaum 38  
22143 Hamburg  
Germany

T: +49 (0)40 675 99 39 0  
F: +49 (0)40 675 99 39 20  
[info@starlab.de](mailto:info@starlab.de)  
[www.starlab.de](http://www.starlab.de)

**STARLAB FRANCE SARL**

30 Rue Jean Rostand  
91400 Orsay  
France

T: +33 (0)1 60 13 71 70  
F: +33 (0)1 69 41 48 65  
[info@starlab.fr](mailto:info@starlab.fr)  
[www.starlab.fr](http://www.starlab.fr)

**STARLAB S.r.l.**

Via Pinturicchio 1  
20133 Milano  
Italy

T: +39 02 7020 1040  
F: +39 02 7020 1033  
[info@starlab.it](mailto:info@starlab.it)  
[www.starlab.it](http://www.starlab.it)